

Press release
No. 031/2020

New countertop disinfectors from Miele for surgeries, chiropodists and tattoo studios

- ▶ Compact, flexible and ergonomic
- ▶ Safe reprocessing combined with reduced water and energy consumption

Gütersloh, June 15, 2020 – There is not always sufficient space for the machine-based reprocessing of instruments, ensuring reliable results using standardised programmes, temperatures and holding times, where disinfected supplies are required. To address this issue, Miele is now offering a new and compact washer-disinfector specifically designed for ENT, dermatology, gynaecology and chiropody surgeries as well as tattoo studios, which sits neatly on a worktop.

Top-notch capacity on the smallest of footprints: The new PWD 8534 model can be fitted with two different baskets and various inserts to suit individual needs. These accommodate instruments from medical and chiropody practices as well as the equipment used by tattoo studios for proper reprocessing. The capacity per cycle, for instance for ENT instruments: Up to 8 lumened instruments as well as further instruments in various inserts.

The new model is convenient to use. Tailored programmes for specific professions can be selected using direct-access buttons. These include for example specific cycles for gynaecological or surgical instruments. The simple logics of the touch-on-glass colour display simplifies the user interface. Visual control over reprocessing is facilitated by a full glass door which, like the glass fascia with its colour display, is easy to clean. Automatic dispensing is via two pumps which deliver liquid process chemicals as and when needed in the relevant phases of the cycle from integrated supply containers.

A built-in water softener ensures top-class cleaning results. Triple filtration is another key contributing factor. The durability of these machines is guaranteed not least by the use of the high-quality stainless steel from which the spray arms, filters and chamber are made. The latter offers a useable volume of 50 l and hence ample space for all types of medical instruments. The DryPlus active hot-air drying of lumens protects against corrosion. A highly efficient HEPA Class H14 particulate filter prevents air-borne particles from ambient room air from entering the drying system and contaminating instruments.

Compliance with process parameters (in particular pressure and temperature) is monitored by two independent and redundant systems. All reprocessing data can either be printed out

Free for
editorial use

Please submit
copies to
Miele & Cie. KG
Public Relations

Postal address
P.O.B.
33325 Gütersloh
Germany

Phone
+49 5241 89-1953
Fax
+49 5241 89-1950

www.miele.de
presse@miele.de

or collected and archived using documentation software in compliance with all legal requirements.

Media contact:

Anke Schläger

Phone: +49 5241 89-1949

Email: anke.schlaeger@miele.com

Company profile: Miele is the world's leading manufacturer of premium domestic appliances including cooking, baking and steam-cooking appliances, refrigeration products, coffee makers, dishwashers and laundry and floor care products. This line-up is augmented by dishwashers, washing machines and tumble dryers for commercial use as well as washer-disinfectors and sterilisers for use in medical and laboratory applications (Miele Professional). The Miele company, founded in 1899, has eight production plants in Germany as well as one plant each in Austria, the Czech Republic, China, Romania and Poland. These are joined by the two plants belonging to Miele's Italian medical technology subsidiary Steelco. Turnover in the 2018/19 business year was around € 4.16 bn, with a share of business outside Germany of approx. 71%. Miele is represented with its own sales subsidiaries and via importers in almost 100 countries. Throughout the world, the family-run enterprise, now in its fourth generation, employs a workforce of around 20,500, whereof 11,050 are employed in Germany. The company has its headquarters in Gütersloh in Westphalia.

There is one photograph with this text



Photo 1: A compact washer-disinfector from Miele which sits on a countertop is now available to meet the requirements of medical surgeries, chiropodists and even tattoo studios. (Photo: Miele)

Text and photo download: www.miele-presse.de

Follow us on:

 @Miele_Press

 @Miele

 @Miele_com

 Miele | Miele Professional